


Fluent™

TRANSITIONLESS
MICRO-INTRODUCER

UP TO 57% LOWER
initial force than the leading competitor.*



Galt Gold™



Experience a **Fluent™** Introduction

Galt Medical Corp., one of the premier providers of .018" micro-access products, is pleased to introduce the Fluent™ Transitionless Micro-Introducer (TMI). The Fluent™ TMI is designed for ease of access, with redesigned angles and silicone coating. Test data shows the Fluent™ is more effective in reducing penetration forces than the leading competitor.

- **Fluent transition**
- **4Fr and 5Fr sizes**
- **Regular and Stiff cannula configurations**
- **X-Sharp™ 21ga echogenic needle with bevel indicator**
- **Variety of guidewire options including high-radiopacity tips**
- **All-inclusive packaging**



*1:1 simulated use test data on file at Galt Medical Corp.

For more information about the Fluent Transitionless Micro-Introducer or to set up a meeting with a Galt Sales Consultant call (800) 639-2800

Fluent™

TRANSITIONLESS MICRO-INTRODUCER

Fluent™ Transitionless Micro-Introducer Kit

| Catalog Number | Size | Description | Guidewire Length | Guidewire Type | Needle |
|----------------|------|-------------|------------------|--------------------|--------|
| KIT-038-01 | 4F | 10cm | 40cm | SS/Galt Gold™ | Echo |
| KIT-038-02 | 5F | 10cm | 40cm | SS/Galt Gold™ | Echo |
| KIT-038-03 | 4F | 10cm | 40cm | Nitinol/Galt Gold™ | Echo |
| KIT-038-04 | 5F | 10cm | 40cm | Nitinol/Galt Gold™ | Echo |
| KIT-038-05 | 4F | 10cm | 60cm | Nitinol/Galt Gold™ | Echo |
| KIT-038-06 | 5F | 10cm | 60cm | Nitinol/Galt Gold™ | Echo |
| KIT-038-07 | 4F | 10cm | 40cm | SS/SS | Echo |
| KIT-038-08 | 5F | 10cm | 40cm | SS/SS | Echo |
| KIT-038-11 | 4F | 10cm | 40cm | Nitinol/SS | Echo |
| KIT-038-12 | 5F | 10cm | 40cm | Nitinol/SS | Echo |
| KIT-038-13 | 4F | 10cm | 40cm | Nitinol/Pt | Echo |
| KIT-038-14 | 5F | 10cm | 40cm | Nitinol/Pt | Echo |
| KIT-038-15 | 4F | 10cm | 40cm | SS/Pt | Echo |
| KIT-038-16 | 5F | 10cm | 40cm | SS/Pt | Echo |

Stiffen Fluent™ Transitionless Micro-Introducer Kit

| Catalog Number | Size | Description | Guidewire Length | Guidewire Type | Needle |
|----------------|------|-------------|------------------|--------------------|--------|
| KIT-039-01 | 4F | 10cm | 40cm | SS/Galt Gold™ | Echo |
| KIT-039-02 | 5F | 10cm | 40cm | SS/Galt Gold™ | Echo |
| KIT-039-03 | 4F | 10cm | 40cm | Nitinol/Galt Gold™ | Echo |
| KIT-039-04 | 5F | 10cm | 40cm | Nitinol/Galt Gold™ | Echo |
| KIT-039-05 | 4F | 10cm | 40cm | SS/SS | Echo |
| KIT-039-06 | 5F | 10cm | 40cm | SS/SS | Echo |
| KIT-039-07 | 4F | 10cm | 40cm | Nitinol/SS | Echo |
| KIT-039-08 | 5F | 10cm | 40cm | Nitinol/SS | Echo |
| KIT-039-11 | 4F | 10cm | 40cm | Nitinol/Pt | Echo |
| KIT-039-12 | 5F | 10cm | 40cm | Nitinol/Pt | Echo |
| KIT-039-13 | 4F | 10cm | 60cm | Nitinol/Galt Gold™ | Echo |
| KIT-039-14 | 5F | 10cm | 60cm | Nitinol/Galt Gold™ | Echo |



Galt Medical Corp. • Garland, Texas • Phone: (800) 639-2800 • Fax (972) 271-4706 • www.galtmedical.com

ISO 13485 Certified. MKT-100-07 Rev. A.

For more information about the Fluent Transitionless Micro-Introducer or to set up a meeting with a Galt Sales Consultant call (800) 639-2800