


# GALTSTICK™ INTRODUCER SYSTEM



Now available with  
**Lubricity+**  
HYDROPHILIC COATING



Radiopaque marker band  
on all sheaths



Kits available  
with guidewire

**The GaltStick™ Introducer System** is a 6 French Coaxial sheath assembly with stiffening stylet. The GaltStick™ allows for the introduction of a .018" guidewire and .038" working guidewire in percutaneous procedures.

- Embedded radiopaque marker band on all sheaths for enhanced visibility
- Smooth dilator/sheath transition
- Variety of kit offerings
- Hydrophilic or silicone coated sheath configurations
- 6Fr outer sheath, 4Fr dilator, 20cm usable length



For more information about the GaltStick or to set up a meeting with a Galt Sales Consultant call (800) 639-2800

## GaltStick™ Introducer System



### GaltStick™ Introducer Set with Radiopaque Band

Model Number	Description	Quantity
KIT-030-04	GaltStick™ Introducer Set with (1) .018" x 60cm Stainless Steel Mandrel/Platinum Coil Guidewire (1) .038" x 145cm PTFE Double Flex Tips/3mm "J" - Straight	1
KIT-030-05	GaltStick™ Introducer Set with (1) .018" x 60cm Nitinol Mandrel/Platinum Coil Guidewire (1) .038" x 145cm PTFE Double Flex Tips/3mm "J" - Straight	1
KIT-030-06	GaltStick™ Introducer Set with (1) .018" x 60cm Stainless Steel Mandrel/Platinum Coil Guidewire	1
INT-021-06	Galt Stick™ Introducer Set Only	1

### Hydrophilic Coated GaltStick™ Introducer Set with Radiopaque Band

Model Number	Description	Quantity
KIT-036-04	GaltStick™ Introducer Set with (1) .018" x 60cm Stainless Steel Mandrel/Platinum Coil Guidewire (1) .038" x 145cm PTFE Double Flex Tips/3mm "J" - Straight	1
KIT-036-05	GaltStick™ Introducer Set with (1) .018" x 60cm Nitinol Mandrel/Platinum Coil Guidewire (1) .038" x 145cm PTFE Double Flex Tips/3mm "J" - Straight	1
KIT-036-06	GaltStick™ Introducer Set with (1) .018" x 60cm Stainless Steel Mandrel/Platinum Coil Guidewire	1
INT-034-06	GaltStick™ Introducer Set Only	1

### Accessories

Model Number	Description	Quantity
NDL-101-00	21ga x 15cm Trocar Needle	5
SGW-015-15	.018" x 60cm Stainless Steel Mandrel Wire - Platinum Coil Guidewire	5
SGW-019-12	.018" x 60cm Nitinol Mandrel Wire - Platinum Coil Guidewire	5

GaltStick™ Specifications: 6 French Sheath, 4 French Dilator, 20cm length

