


MICRO-ACCESS/MST INTRODUCER KITS AND SETS

Complete MST Kit



Smooth transition



Tearaway introducer available in half-sizes



Kink-resistant core-wire design



Includes safety needle



Introducer Sets available

GALT's Micro-Access/MST Introducer Kits and Sets were precisely designed for smaller access sites. Featuring 5cm, 7cm, and 10cm lengths with the option of half-sized configurations. Each introducer features a smooth transition, easy and secure locking system, unique wing design, and funnel-shaped lead-in.

- Full kits and introducer sheath only
- 5cm, 7cm, and 10cm lengths
- Half-sized configurations
- All kits include a Nitinol guidewire
- 10cm lengths include marked Nitinol wire
- Safety components included



Micro-Access Tearaway Introducer Kits

Kit contents: PTFE Tearaway Introducer Set, a Nitinol Mandrel Guidewire/SS Coil, a safety guidewire introducer needle, a #11 safety scalpel, and a standard guidewire introducer needle

Kits contain a .018" x 40cm Marked Nitinol Mandrel Wire with a Stainless Steel Soft Tip

Catalog Number	Size	Description	Needle	Units/Box
Kit-010-00	3.0F	5.0cm	21ga. x 5cm. Safety	5
Kit-010-01	3.5F	5.0cm	21ga. x 5cm. Safety	5
Kit-010-02	4.0F	5.0cm	21ga. x 5cm. Safety	5
Kit-010-03	4.5F	5.0cm	21ga. x 5cm. Safety	5
Kit-010-04	5.0F	5.0cm	21ga. x 5cm. Safety	5
Kit-010-05	5.5F	5.0cm	21ga. x 5cm. Safety	5
Kit-010-06	6.0F	5.0cm	21ga. x 5cm. Safety	5

Kits contain a .018" x 65cm Marked Nitinol Mandrel Wire with a Stainless Steel Soft Tip

Catalog Number	Size	Description	Needle	Units/Box
Kit-010-10	3.0F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-11	3.5F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-12	4.0F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-13	4.5F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-14	5.0F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-15	5.5F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-16	6.0F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-17	7.0F	10.0cm	21ga. x 7cm. Safety	5
Kit-010-18	8.0F	10.0cm	21ga. x 7cm. Safety	5

Micro-Access Tearaway Introducer Sets

Set contents: PTFE Sheath and Dilator only

5cm Introducer Sets

Order Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-106-00	3.0F	5cm	.018"	10
INT-106-01	3.5F	5cm	.018"	10
INT-106-02	4.0F	5cm	.018"	10
INT-106-03	4.5F	5cm	.018"	10
INT-106-04	5.0F	5cm	.018"	10
INT-106-05	5.5F	5cm	.018"	10
INT-106-06	6.0F	5cm	.018"	10

7cm Introducer Sets

Order Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-107-30	3.0F	7cm	.018"	10
INT-107-40	4.0F	7cm	.018"	10
INT-107-50	5.0F	7cm	.018"	10
INT-107-60	6.0F	7cm	.018"	10
INT-107-70	7.0F	7cm	.018"	10

10cm Introducer Sets

Order Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-106-10	3.0F	10cm	.018"	10
INT-106-11	3.5F	10cm	.018"	10
INT-106-12	4.0F	10cm	.018"	10
INT-106-13	4.5F	10cm	.018"	10
INT-106-14	5.0F	10cm	.018"	10
INT-106-15	5.5F	10cm	.018"	10
INT-106-16	6.0F	10cm	.018"	10
INT-106-17	7.0F	10cm	.018"	10
INT-106-18	8.0F	10cm	.018"	10
INT-106-20	6.5F	10cm	.018"	10

Long Introducer Sets

Order Number	Introducer Size	Introducer Length	Accepts Wire Diameter	Units/Box
INT-106-31	5.0F	30cm	.018"	5
INT-106-33	5.0F	50cm	.018"	5
INT-106-37	7.0F	50cm	.018"	5

