


MICRO-INTRODUCER KITS & STIFFEN MICRO-INTRODUCER KITS





Galt's Micro-Introducer Kits provide the high quality and precision design expected in a fine medical device. Each feature was designed to aid in the access process, allowing for consistent performance every time.

Galt's line of Stiffen Micro-Introducers sets a new standard for high quality and precision design in micro-introducers. Kits combine our precise dilator and sheath transitions with a reinforced assembly to facilitate smooth entry into difficult vessels.

4 and 5 French sizes available
Variety of wire configurations

Standard or Echogenic needles
All-inclusive packaging



Micro-Introducer Kits

.018" Diameter Stainless Mandrel - Stainless Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 7cm
KIT-002-00	4.0F	10cm	40cm	Standard
KIT-002-01	5.0F	10cm	40cm	Standard
KIT-002-30	4.0F	10cm	40cm	Echogenic
KIT-002-31	5.0F	10cm	40cm	Echogenic
KIT-002-36	4.0F	10cm	60cm	Echogenic
KIT-002-37	5.0F	10cm	60cm	Echogenic

.018" Diameter Stainless Mandrel - Platinum Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 7cm
KIT-002-28	4.0F	10cm	40cm	Standard
KIT-002-29	5.0F	10cm	40cm	Standard
KIT-002-56	4.0F	10cm	60cm	Standard
KIT-002-57	5.0F	10cm	60cm	Standard

.018" Diameter Nitinol Mandrel - Stainless Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 7cm
KIT-002-08	4.0F	10cm	40cm	Standard
KIT-002-09	5.0F	10cm	40cm	Standard
KIT-002-32	4.0F	10cm	40cm	Echogenic
KIT-002-33	5.0F	10cm	40cm	Echogenic
KIT-002-38	4.0F	10cm	60cm	Echogenic
KIT-002-39	5.0F	10cm	60cm	Echogenic

.018" Diameter Nitinol Mandrel - Platinum Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 7cm
KIT-002-12	4.0F	10cm	40cm	Standard
KIT-002-13	5.0F	10cm	40cm	Standard
KIT-002-34	4.0F	10cm	40cm	Echogenic
KIT-002-35	5.0F	10cm	40cm	Echogenic
KIT-002-40	4.0F	10cm	60cm	Echogenic
KIT-002-41	5.0F	10cm	60cm	Echogenic

Micro-Introducer Kits with Side Holes

.018" x 40cm Nitinol Mandrel - Stainless Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 7cm
KIT-033-04	4.0F	10cm	40cm	Echogenic

Stiffen Micro-Introducer Kits

.018" Diameter Stainless Mandrel - Stainless Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 12cm
KIT-011-00	4.0F	10cm	40cm	Standard
KIT-011-01	5.0F	10cm	40cm	Standard
KIT-011-02	4.0F	10cm	40cm	Echogenic
KIT-011-03	5.0F	10cm	40cm	Echogenic

.018" Diameter Stainless Mandrel - Platinum Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 12cm
KIT-011-20	4.0F	10cm	40cm	Standard
KIT-011-22	4.0F	10cm	40cm	Echogenic
KIT-011-23	5.0F	10cm	40cm	Echogenic

.018" Diameter Nitinol Mandrel - Stainless Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 12cm
KIT-011-40	4.0F	10cm	40cm	Standard
KIT-011-41	5.0F	10cm	40cm	Standard
KIT-011-42	4.0F	10cm	40cm	Echogenic
KIT-011-43	5.0F	10cm	40cm	Echogenic

.018" Diameter Nitinol Mandrel - Platinum Tip Wire

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 12cm
KIT-011-60	4.0F	10cm	40cm	Standard
KIT-011-61	5.0F	10cm	40cm	Standard
KIT-011-62	4.0F	10cm	40cm	Echogenic
KIT-011-63	5.0F	10cm	40cm	Echogenic
KIT-011-66	4.0F	15cm	45cm	Standard
KIT-011-67	5.0F	15cm	45cm	Standard

Bariatric Micro-Introducer Kits

.018" x 60cm Nitinol Mandrel - Platinum Tip Wire*

Order Number	Introducer Size	Introducer Length	Wire Length	Needle 21ga x 12cm
KIT-011-99	4.0F	15cm	60cm	Standard

*Note: KIT-011-99 contains a 19cm dilator and 12cm needle.



NOTE: Order Quantity - 10 kits/box