

TEARAWAY INTRODUCER SETS & TEARAWAY INTRODUCER KITS



Unique wing design



Easy-to-thread funnel system



Patented locking system



Needle bevel indicator



Smooth transition

GALT's Tearaway Introducer Sets & Kits combine precision ergonomic design with the high quality expected in fine medical devices. Each introducer system provides consistent performance and offers a variety of kit configurations for every need. The smooth transition, easy and secure lock, funnel lead-in and unique patented ergonomic hub design add to a successful access every time.

- 14cm, 30cm, 44cm, and 50cm lengths available
- Sizes available from 4Fr-16Fr, including some half-sizes
- Full kits or individual sets



Tearaway Introducer Kits

Kit contents: a 14cm Tearaway Introducer Set, 50cm flex tip/3mm "J"-Straight Stainless Guidewire, a 7cm introducer needle, and a 10cc Luer slip syringe

14cm Tearaway Introducer Kits

Catalog Number	Introducer Size	Accepts Wire Diameter	Length	Units/Box
KIT-009-01	4.0F	.028"	50cm	5
KIT-009-02	5.0F	.038"	50cm	5
KIT-009-03	6.0F	.038"	50cm	5
KIT-009-04	7.0F	.038"	50cm	5
KIT-009-05	8.0F	.038"	50cm	5
KIT-009-07	9.0F	.038"	50cm	5
KIT-009-08	9.5F	.038"	50cm	5
KIT-009-09	10.0F	.038"	50cm	5
KIT-009-10	10.5F	.038"	50cm	5
KIT-009-11	11.0F	.038"	50cm	5
KIT-009-12	12.0F	.038"	50cm	5
KIT-009-13	12.5F	.038"	50cm	5
KIT-009-14	13.0F	.038"	50cm	5
KIT-009-15	14.0F	.038"	50cm	5
KIT-009-16	15.0F	.038"	50cm	5
KIT-009-17	16.0F	.038"	50cm	5

5cm, 7cm, and 10cm Micro-Access Introducer Kits also available

Tearaway Introducer Sets

Set contents: PTFE Sheath and Dilator only

14cm Tearaway Introducer Sets

Catalog Number	Introducer Size	Accepts Wire Diameter	Units/Box
INT-105-02	4.0F	.028"	10
INT-105-04	5.0F	.038"	10
INT-105-06	6.0F	.038"	10
INT-105-07	7.0F	.038"	10
INT-105-08	8.0F	.038"	10
INT-105-09	9.0F	.038"	10
INT-105-10	9.5F	.038"	10
INT-105-11	10.0F	.038"	10
INT-105-12	10.5F	.038"	10
INT-105-13	11.0F	.038"	10
INT-105-14	12.0F	.038"	10
INT-105-15	12.5F	.038"	10
INT-105-16	13.0F	.038"	10
INT-105-17	14.0F	.038"	10
INT-105-18	15.0F	.038"	10
INT-105-19	16.0F	.038"	10

Long Length Tearaway Introducer Sets

Catalog Number	Introducer Size	Accepts Wire Diameter	Length	Units/Box
INT-105-21	5.0F	.038"	30cm	5
INT-105-25	7.0F	.038"	30cm	5
INT-105-29	8.0F	.038"	30cm	5
INT-105-33	9.0F	.038"	30cm	5
INT-105-34	9.0F	.038"	44cm	5
INT-105-37	10.0F	.038"	30cm	5
INT-105-41	11.0F	.038"	30cm	5
INT-105-45	12.0F	.038"	30cm	5
INT-105-46	12.0F	.038"	44cm	5
INT-105-49	14.0F	.038"	30cm	5
INT-105-53	16.0F	.038"	30cm	5
INT-105-55	16.0F	.038"	50cm	5

5cm, 7cm, and 10cm Micro-Access Introducer Sets also available

