



MicroSlide™

Pediatric MST Products

For smaller patient access, keep it short and sweet.

Galt Medical Corp. introduces the MicroSlide™ line of vascular access products for smaller patients and smaller vessels. Maintaining the features that are preferred in our current offerings, this gives the clinician an added advantage to size the product to the patient. Patients benefit from being accessed by the appropriately sized products. It is a win-win solution, and now available from Galt Medical.

2Fr x 3cm and 3Fr x 3cm tearaway introducers

Available individually or in complete MST Kit

.008" and .010" diameter guidewires

24ga x 2.5cm and 24ga x 4cm needles



(Actual Size)




For more information on the features and benefits of Galt Medical Corp. products, please visit www.GaltMedical.com or call 800.639.2800

MicroSlide™

Pediatric MST Products

Growth Stages

The MicroSlide™ Pediatric MST product line is the first to offer both a 2Fr and 3Fr MST offering. Galt continues to offer the flexibility of a full MST Kit or sterile, separate components. Our MST Tearaway line offering continues to grow and now includes 3cm, 5cm, 7cm and 10cm lengths. Our popular .008" and .010" guideline configurations are now offered in 20cm, 40cm and 60cm lengths. We are committed to continue this growth with a priority on patient care.



(Actual Size)

Tearaway Introducer Sets

Set Contents: PTFE Sheath and Dilator only

Catalog Number	French Size	Length	Accepts Guidewire	Units/Box
INT-101-02	2F	3cm	.008"	10
INT-101-03	3F	3cm	.010"	10



Tearaway Introducer Kits

Catalog Number	French Size	Length	Guidewire	Needle	Units/Box
KIT-051-00	2F	3cm	.008" x 20cm	24ga x 2.5cm	5
KIT-051-02	2F	3cm	.008" x 20cm	24ga x 4cm	5
KIT-051-11	3F	3cm	.010" x 40cm	24ga x 4cm	5



Guidewires

All guidewires are nitinol mandrel-tapered both ends

Catalog Number	Size	Length	Units/Box
SGW-051-06	.008"	20cm	10
SGW-051-08	.008"	40cm	10
SGW-051-03	.008"	60cm	10
SGW-051-07	.010"	20cm	10
SGW-051-09	.010"	40cm	10
SGW-051-04	.010"	60cm	10



(Actual Size)

Needles

Catalog Number	Size	Length	Tip	Units/Box
NDL-100-29	24ga	2.5cm	echogenic	25
NDL-100-30	24ga	4cm	echogenic	25

